

US Offices

Corporate Headquarters

Caliper Life Sciences, Inc.
68 Elm Street
Hopkinton, MA 01748-1668
Telephone: 508-435-9500
Fax No: 508-435-3439
Website: www.caliperLS.com

Worldwide Offices

France

Caliper Life Sciences S.A.
ZAC PARIS-NORD II
13 rue de la Perdrix
BP 40016
95911 Roissy Charles de Gaulle Cedex
France
Telephone: +33-1-48-63-71-35
Fax No: +33-1-48-63-71-53

Benelux

Caliper Life Sciences NV/SA
Klapstraat 13
B-1790 Teralfene
Belgium
Telephone: +32-53-66-26-70
Fax: +32-53-66-27-32

Germany

Caliper Life Sciences GmbH
Eisenstr. 9c
DE-65428 Rüsselsheim
Germany
tel. +49 6142 834 93-0
fax: +49 6142 162 821

Japan

Caliper Life Sciences Japan
Saito-Bldg. 2F
Yushima 2-17-15, Bunkyo-ku
Tokyo 113-0034
Japan
Telephone: + 81-3-5840-6551
Fax: + 81-3-5840-6554

Switzerland

Caliper Life Sciences AG
Nordstrasse 17
CH-4665 Oftringen
Switzerland
Telephone: + 41-62-788-7000
Fax: + 41-62-788-7017

United Kingdom

Caliper Life Sciences Ltd.
1 Wellfield
Preston Brook
Runcorn, Cheshire WA7 3AZ
United Kingdom
Telephone: +44-1928-711448
Fax No: +44-1928-791228

©2004 Caliper Life Sciences, Inc.
Caliper and LabChip are registered trademarks, and the Caliper Logo and the phrase "Working Innovation" are trademarks of Caliper Life Sciences, Inc. Zymark and AutoTrace are registered trademarks of Zymark Corporation.

AP23 03/04

Conclusion

The AutoTrace is specifically designed for extraction of large volume water samples. The extraction of silica gel treated hexane extractable material (SGT-HEM) is efficient using Isolute SPE column TPH. The positive flow feature of the AutoTrace ensures that each extraction is uniform and optimizes EPA Method 1664.

Extraction of Total Petroleum Hydrocarbons Using the AutoTrace

IST International Sorbent Technology, Dyffryn Business Park, Hengoed, Mid Glamorgan, CF82 7RJ, United Kingdom

Lynn Jordan and Ed Alderman

Caliper Life Sciences Corporation, Applications Laboratory, Hopkinton, Massachusetts, USA

Overview

This method was developed for the extraction of total petroleum hydrocarbons (or silica gel treated-hexane extractable material (SGT-HEM) as described in EPA method 1664) from water samples. The analytes are determined by gravimetric analysis.

Extraction Procedure

*ISOLUTE® SPE Column:	TPH 1 g / 6 ml
Pre-treatment:	Add 10 mL of methanol to 1 L of sample. Acidify the sample with 6 M HCl to a pH ~2.
Solvation:	Program the AutoTrace to WASH SYRINGE with 5 mL of methanol, followed by a CONDITION COLUMN with 10 mL of methanol.
Equilibration:	Program the AutoTrace to WASH SYRINGE with 10 mL of distilled deionised water pH~2, followed by a CONDITION COLUMN with 10 mL of distilled deionised water pH~2.
Sample application:	Program AutoTrace to LOAD 1100 mL SAMPLE. Following the sample load step, include a PAUSE AND ALERT step.
Interference elution:	During pause step, remove sample line from sample bottle. Manually rinse sample bottle with 10 mL of acetone, swirl well to cover sides of bottle. Dilute acetone in bottle with 40 mL of distilled deionised water, pH~2. Return sample line to bottle, and press CONT on AutoTrace. Program AutoTrace to LOAD 60 mL SAMPLE. Following this step, include a PAUSE AND ALERT step. Repeat this step if necessary until the rinse solvent appears clear. Press CONT on AutoTrace. Dry cartridge assembly for 30 min.
Analyte elution:	Rinse AutoTrace syringe with 5 mL of hexane. ELUTE TO SOAK with 4 mL of hexane. Follow this step with a TIMED PAUSE for 2 min. ELUTE TO COLLECT with 4 mL of hexane into tared collection tubes. Concentrate solvent to near dryness in a gentle stream of nitrogen. Collection tubes may be placed in a heating block held at 35 C to expediate evaporation. Once solvent is almost gone, weigh tubes in one minute intervals until weight loss is less than 1 mg. Total petroleum hydrocarbons is the weight of residue.

Setup parameters on AutoTrace should be as follows:

Cond Flow:	40.0 mL / min	Push Delay	5 sec
Load Flow:	30.0 mL / min	Air Factor	1.0
Rinse Flow:	40.0 mL / min	Autowash Vol	1.00 mL
Elute Flow:	20.0 mL / min		
Cond Air Push:	40.0 mL / min		
Rinse Air Push:	40.0 mL / min		
Elute Air Push:	40.0 mL / min		

* ISOLUTE column part numbers represent the product configuration of choice for use with an AutoTrace SPE Workstation.

The AutoTrace procedure should be written as follows:

- Step 1: Process 6 samples using the following procedure:
- Step 2: Wash syringe with 5 mL of methanol
- Step 3: Condition column with 10 mL of methanol in SOLVENT WASTE
- Step 4: Wash syringe with 10 mL of water pH=2
- Step 5: Condition column with 10 mL of water to AQUEOUS WASTE
- Step 6: Load 1100 mL of sample on column
- Step 7: Pause and Alert operator, resume when CONTinue is pressed
- Step 8: Load 60 mL of sample onto column
- Step 9: Pause and Alert operator, resume when CONTinue is pressed
- Step 10: Dry column with gas for 30 minutes
- Step 11: Wash syringe with 5 mL of hexane
- Step 12: Soak and Collect 4 mL fraction using hexane
- Step 13: Pause for 2 minutes
- Step 14: Collect 3 mL fraction into sample tube using hexane
- Step 15: Clean each sample path with 10 mL into SOLVENT WASTE
- Step 16: Clean each sample path with 50 mL into AQUEOUS WASTE
- Step 17: Wash syringe with 10 mL of water
- Step 18: END

NOTE: Place sample lines into methanol for step 15 and reagent water for step 16 to flush out lines.

Structural Various	Non-polar hydrocarbons.
Structural Considerations	This method is suitable for the non-polar petroleum hydrocarbons (SGT-HEM).
Matrix Considerations	The matrix is polar, and the analytes are extracted by a non-polar retention mechanism.
Analytical Method	Gravimetric analysis using an analytical balance having precision to 0.1 mg.
General Comments	<ol style="list-style-type: none"> 1. This method describes an automated procedure for the determination of silica gel treated-hexane extractable material (SGT-HEM) from an aqueous sample. This column can also be used for Total Oil and Grease type extractions (EPA 1664 HEM) if the concentration of polar components is low. In this situation, the the hexane elution step should be followed by an elution (into a second tared vial) using 2 x 4 mL THF/hexane (1:1, v/v).Combination of these two fractions will give a determination of Total Oil and Grease, as described in EPA method 1664.However, for higher concentrations of polar components, we would recommend reference to application note IST1018 (TPH/Total Oil and Grease determination on the AutoTrace) which utilises a column specifically optimised for EPA 1664 extractions on the AutoTrace system. 2. Due to the nature of the analytes, the bottle washing steps after sample loading are very important, as analytes do stick to the walls of the sample bottle. For this reason, sample splitting is not recommended for TPH or Oil and Grease samples.