

## Minimizing fluorescence using a 457nm laser excitation wavelength

### Introduction <Benefits of Raman Spectroscopy>

Raman spectroscopy is a popular method for analyzing molecular structure and is considered a complementary technique to infrared spectroscopy. Recently, Raman has been attracting the attention of FT-IR users because it offers several important advantages over FT-IR. Raman spectroscopy is a non-contact and non-destructive technique and measurements can be made with little to no sample preparation. Samples can be measured with a spatial resolution as small as 1 $\mu$ m and depth profiling can also be easily performed on transparent samples. However, in some cases, good quality analysis by Raman spectroscopy can be adversely affected by interference from fluorescence.

The JASCO NRS series of Raman spectrometers incorporates several fluorescence compensation features, including; a confocal optical system, multiple laser excitation wavelengths and a simple, but powerful, fluorescence elimination algorithm (patented).

In this application we evaluated the NRS-4100 (Figure 1) fitted with a 457 nm laser for the measurement of samples that exhibit fluorescence, to determine if this can be a better alternative to red and NIR lasers such as 785 or 1064 nm excitation.



Fig. 1 Raman Spectrometer

### Methods for minimizing fluorescence

Both Raman scattering and fluorescence are phenomena where the wavelengths of light emitted from a sample are different than the wavelength of the input excitation light. If the wavelengths of Raman scattering and fluorescence overlap, it is impossible to obtain good Raman spectra. The NRS series Raman instruments have three tools for reducing or eliminating unwanted fluorescence.

1. An important method for minimizing fluorescence is the fluorescence correction algorithm (JASCO patent) included with the Spectra Manager II software.
2. When fluorescence is emitted from the surrounding matrix, the most widely used method to reduce this fluorescence is to use higher spatial resolution by reducing the aperture in the confocal optical system.
3. When the sample itself emits fluorescence, the most effective method is to change the wavelength of the excitation laser. The wavelengths of Raman scattering do not change even if the excitation wavelength is changed, while the fluorescence wavelengths are dependent on the excitation wavelength. Therefore, it is possible to minimize or even eliminate the overlap of the Raman scattering and fluorescence by changing the excitation wavelength. As Raman spectra are displayed as a shift value from the excitation wavelength, Raman peaks always appear at the same position independent of the excitation wavelength and so Raman spectra with substantially minimized fluorescence can be obtained. An additional benefit of using the 457 nm laser is that the intensity of Raman scattering is inversely proportional to the fourth power of excitation wavelength, therefore by using a shorter wavelength laser as the excitation source, Raman scattering intensity can be substantially increased while minimizing fluorescence.




Fig. 2 457 nm laser spot  
(a blue-purple-color)

The 457 nm laser as shown in Figure 2 has a shorter wavelength than the 532 nm laser and offers two important advantages; the same CCD detector can be used as with the 532 nm laser and the Raman scattering intensity can be up to 1.8 times higher when using a laser of equal power output.

## System configuration

- NRS-4100 Raman spectrometer
- 532 nm laser (100 mW)
- 785 nm laser (100 mW)
- 457 nm laser
- Laser switching mechanism
- 900 and 400 gr/mm grating

## Measurement and Analysis

A fiber sample was measured using the laser excitation wavelengths of 457 nm, 532 nm and 785 nm; measurement parameters are outlined in table 1. The spectra in Figure 3 demonstrate that using excitation wavelengths at both 532 nm and 785 nm result in strong fluorescence which completely obscures any peaks of the Raman spectrum. However, by using an excitation wavelength at 457 nm, the fluorescence interference exhibited by the other excitation lasers is quite reduced. The spectrum obtained using the 457 nm laser and the fluorescence correction software resulted in data which could be compared to a database library spectrum and the measured sample was correctly identified as a nylon-6 fiber (Figure 4).

Table 1: Excitation Wavelength Measurement Parameters

Wavelength [nm]	Grating [Line/mm]	Measurement time[sec]	Number of scans
457	900	30	2
532	900	5	12
785	400	30	2




Fig. 3 Overlaid spectra for three laser excitation wavelengths




Fig. 4 Database Search Result



# Application Note

260-AN-0012

---

## Summary

A nylon sample fiber was measured using the 'standard' 532 and 785 nm excitation wavelengths and the results for both lasers demonstrate a strong fluorescence emission. However, using the shorter (higher energy) 457 nm laser excitation, which offers a significantly higher Raman scattering intensity, proves to be an extremely effective method to minimize the effects of fluorescence.

In addition to the nylon test fiber, we used the three-laser NRS-4100 Raman spectrometer (582 nm, 457 nm and 785 nm) to evaluate a range of samples that exhibit strong fluorescence, such as polyimide and biological materials and found that the 457 nm excitation wavelength offers spectra with much lower fluorescence than the 'standard' combination of 532 nm and 785 nm laser excitation wavelengths.