

## Measurement of complex forming reaction of Nickel Sulfate and Rochelle salt

Stopped-flow CD measurement technique is well known method to analyze unfolding and refolding process of Protein and also observe complex forming reaction. Because Transition metal complex have typically absorbance in region from Visible to N-IR range.

We introduce this measurement about complex forming reaction of Nickel Sulfate and Rochelle Salt in using of high speed stopped flow systems consists of J-820 and SFS-492.

### <CD spectra of mix solution including Nickel Sulfate and Rochelle Sault>

Fig. 1 shows that CD spectrum of the mixed solution sample which consists of 0.24M Nickel Sulfate solution and 0.36M Rochelle Salt solution as 1:1 mixture ratio. It can show CD signal in broad range from UV/Vis to N-IR.




Fig. 1: CD spectra of mix solution including Nickel Sulfate and Rochelle Salt  
Cell path :10 mm Instrument range: 250-900 nm (J-820) 900-1350 nm (J-730)

### <Stopped- Flow measurement>

The complex forming process about the above sample is measured with CD stopped flow system. Following Fig.2 shows the CD spectra in NIR-region 720nm and Fig.3 shows the one in N-IR region 1000 nm.

It shows in Fig.4 that complex forming reaction is finished within 100msec and shapes of both data in 720 nm, 1000 nm is match after data normalization, which mean both data in 720 nm and 1000 nm indicates same reaction process.




Fig. 3: Complex forming reaction of Nickel Sulfate and Rochelle salt in 1000 nm.




Fig. 4: Stopped flow data comparison between 720 nm (Blue) and 1000 nm (Green)

### <Measurement Condition>

Syringe 1 : 0.24M Nickel Sulfate Solution  
 Syringe 2 : 0.36M Rochelle Salt Solution  
 Mixture ratio : 100  $\mu$ l : 100  $\mu$ l  
 Total flow rate : 5 ml/sec  
 Cell length : 2 mm  
 Measurement range : 720 nm (SBW5 nm, Data acquisition 75 times)  
                           1000 nm (SBW10nm, Data acquisition 50 times)  
 Data pitch : 1 msec  
 Response : 2 msec

### <Reference>

- (1) Hiroyuki Miyake, Hideki Sugimoto, Hitoshi Tamiaki and Hiroshi Tsukube, (2005) *Chem. Commun.*, 4291-4293
- (2) Hiroyuki Miyake, Hiroshi Kamon, Ikuko Miyahara, Hideki Sugimoto, and Hiroshi Tsukube, (2008) *J. Am. Chem. Soc.*, **130**, 792-793