

## Example of the measurement of a highly reflective material, a dielectric multi-layer mirror using absolute reflectance spectroscopy

### Introduction

Dielectric multi-layer mirrors are laminated optics comprising of a combination of high refractive index and low refractive index materials. Using the interference effect, this mirror can provide extremely high reflectance (near 100%) in a specific wavelength range. (Fig. 1)




Fig.1 Diagram of the dielectric multi-layer mirror

This type of mirror is widely used in cameras, telescopes and optics for optical communication to reduce loss in light intensity. Therefore, in the quality control and R&D of dielectric multi-layer mirrors, it is very important to evaluate the reflectance with very high accuracy.

JASCO provides a high performance absolute reflectance measurement system using the V-700 series UV/Vis spectrophotometer with double-beam optical system, which has high photometric stability. The absolute reflectance system can perform measurement at arbitrary and user-selectable incident angles; this allows the measurement of dielectric multi-layer mirrors at designated incident angles.

As an example of high reflectance measurement using the JASCO absolute reflectance measurement system, this application note reports the results for a dielectric multi-layer mirror.

**Keyword:** V-700 UV-visible spectrophotometer, Dielectric multi-layer mirror, Absolute reflectance, High reflectance, Measurement repeatability

### Sample

Dielectric multi-layer mirror

### Measurement

- (1) Dark measurement is performed with a light shielding plate.
- (2) Baseline measurement is performed against air.
- (3) Sample measurement is performed.

The procedure detailed above was repeated three times.

### Measurement Condition

Bandwidth 5.0 nm  
 Scanning speed 20 nm/min  
 Response 3.84 sec  
 Data interval 1 nm  
 Incident angle 10 degree  
 Polarization p-polarized light



Fig.2 Absolute reflectance measurement accessory

## Measurement result

Figure 3 shows an overlay of the absolute reflectance spectra of the dielectric multi-layer film mirror, and Figure 4 shows a zoomed view of the same spectra. As shown in the Figure 4, the measurement repeatability around 100%R is better than 0.1%R, which demonstrates that this system has very high measurement repeatability. Table 1 shows the comparison between the expected values and the measured values. The difference between the average of the measurement value and the expected value is better than 0.15%, which shows that this system has high accuracy and good correlation for reflectance measurement close to 100%R.




Fig. 3 Absolute reflectance spectra of the dielectric multi-layer mirror




Fig. 4 Absolute reflectance spectra of the dielectric multi-layer mirror at around 100%R

Wavelength [nm]	Expected value [%R]	Measurement value [%R]			Average value [%R]	Standard variation	Variation coefficient [%]	Difference between expected and average
		1	2	3				
680	99.9914	99.8663	99.8537	99.8407	99.8536	0.013	0.013	-0.1378
700	99.9918	99.8892	99.8943	99.8578	99.8804	0.020	0.020	-0.1114
720	99.9839	99.9094	99.9316	99.8889	99.9100	0.021	0.021	-0.0739

Table1. Comparison of theoretical and measured values