



## SEROTONIN IN URINE



THE SOUNDEST LC-EC APPLICATIONS FOR  
CLINICAL & DIAGNOSTICS ANALYSIS  
EVER BUILT

Catecholamines  
Serotonin  
Metanephrines  
VMA  
HVA  
5-HIAA  
Homocysteine  
Glutathione  
(di-)sulfides  
Iodide  
Vitamins A, C, D, E, and K  
Q10  
Ubiquinolins



### INTRODUCTION

Serotonin is synthesized by enterochromaffin cells of the intestine and certain neurons of the central nervous system. In blood more than 97% are stored in platelets [1]. Physiological actions of serotonin include the control of circadian rhythms, sleep regulation, sex drive and thermoregulation as well as the influence on melatonin synthesis and on aldosterone regulation [2,3]. Various diseases are related to a pathologic serotonin metabolism [4, 5]. An elevated plasma concentration of serotonin and an increased renal secretion of the serotonin metabolite 5-hydroxyindoleacetic acid may be found in patients with epilepsy. Migraine is associated with a decreased platelet serotonin concentration. Serotonin metabolism is also disturbed in patients who suffer from schizophrenia, autism or psychotic depression. The determination of plasma serotonin level is of decisive importance for the diagnosis of the carcinoid syndrome which is mainly accompanied by an elevated serotonin production [6].

- Standardized, fast and reliable assay
- Kit for standardized sample prep
- Robust & reproducible

### Summary

HPLC with electrochemical detection has been established as a fast and reliable method for the determination of serotonin, catecholamines and metabolites in plasma and urine [12, 13]. The ALEXYS Clinical Analyzer together with a commercially available sample prep kit is dedicated and standardized for routine analysis of urinary serotonin.



Fig. 1. ALEXYS Clinical Analyzer.

## Method

A Recipe ClinRep® complete kit contains all the necessary chemicals and (calibration) materials for sample preparation and analysis. Urine samples are processed as follows:

- 2 mL acidified urine sample (10 mL conc. 32% HCl per liter urine) or urine calibrator is mixed with 4 mL stabilising reagent S and 50 µL internal standard (IS) and subsequently adjusted to a pH 4.5 – 6.5 using 0.5M NaOH.
- The mixture is applied to a ClinRep® sample preparation columns to trap the serotonin present in the sample.
- The column is washed with 15 mL HPLC-grade water and subsequently with 2 mL washing solution W to remove interfering components.
- 5 mL of eluting reagent E is then used to elute serotonin from the extraction column.
- The eluate is collected, mixed (vortex-mixer) and 20 µL injected in the LC system.




Fig. 1 Analysis of 20 µL ClinCal® urine calibrator with a concentration of 346 µg/L serotonin.

The quantification of the serotonin in the urine samples is performed by means of a single-point calibration method using a urine calibrator. The calibrator is a lyophilised urine sample with a known amount of serotonin. The urine calibrator should be processed via the same sample preparation method as the urine samples. An example chromatogram of a urine calibrator analysis is shown in figure 1. An internal standard compensates for recovery losses during the sample preparation step.

Table 1

Set-up	
HPLC	ALEXYS Clinical Analyzer
Flow cell	GC type flow cell with Ag/AgCl saltbridge REF
Column	ClinRep® Analytical column for serotonin in urine

## Analysis of ClinChek® controls

For quality control of the analytical determination Recipe ClinChek® urine controls have been used in both the normal (level I) and the pathological range (level II).




Fig. 2. Overlay of 2 chromatograms of 20 µL injections of ClinChek® control level I (red) and II (blue).

The control samples are lyophilised urine samples which have to be processed in the same way as the urine samples. Both Control I and Control II were analysed and the analyte concentrations quantified using the ClinCal urine calibrator. For both urine controls level I and II the determined serotonin concentrations were within the concentration ranges specified by Recipe on the urine control data sheet (see table I).

Table I. Calculated serotonin concentration in urine controls level I and II,  $n = 4$  (injections)  $\times$  3 (days). Concentration range specified by Recipe is given for reference (source: data sheet supplied with controls).

Component	Specified conc (µg/l)		Calculated conc (µg/l)	RSD (%)
	Min	Max		
Control, level I	123	185	180	1.2
Control, level II	420	630	557	1.0

## Analysis of urine samples

Urine samples of an apparently healthy volunteer were collected and analysed multiple times to determine the recoveries, LOD, intra- and inter-assay precision of the method.

The intra-assay precision of the method was determined using two spiked urine samples A and B, matching the serotonin concentration of control level I and level II, respectively. The urine samples were worked-up 5 times and duplicate analysis were performed to determine the relative standard deviation (RSD, %). This procedure was repeated for 3 days. The RSD's found for sample A and B were smaller than 2%.

Table II. Intra-assay precision for the analysis of serotonin in spiked urine sample A and B,  $n=5$  (samples)  $\times$  2 (injections).

Component	RSD (%)	Conc. ( $\mu\text{g/l}$ )
<i>Sample A</i>		
Day 1	0.6	190
Day 2	1.4	180
Day 3	0.9	177
<i>Sample B</i>		
Day 1	1.1	591
Day 2	1.6	556
Day 3	1.6	551




Fig. 3. Overlay of 10 chromatograms of 20  $\mu\text{L}$  injections of urine sample B.

For all urine samples, controls and calibrator recoveries typically in the range of 80 – 95% were found, compared to a directly injected standard. The concentration limit of detection ( $C_{\text{LOD}}$ ) for

the method was approximately 0.3  $\mu\text{g/L}$  for serotonin. The  $C_{\text{LOD}}$  is calculated based on a 20  $\mu\text{L}$  injection and defined as the concentration that gives a signal that is three times the peak-to-peak noise. The method is linear for the determination of serotonin in the concentration range from 1 – 1000  $\mu\text{g/L}$  [14].

The inter-assay precision of the method was determined over a time period of three days for sample A and B. Both samples were worked-up 5 times and analysed (duplicate injection) every single day and the relative standard deviation calculated.

Table III. Inter-assay precision for the analysis of serotonin in sample A and B.  $n=5$  (samples)  $\times$  2 (duplicate injections)  $\times$  3 (days).

Component	RSD (%)	Conc. ( $\mu\text{g/l}$ )
Sample A	3.2	182
Sample B	3.5	566

The RSD's for the analysis of sample A & B were smaller than 4%.



## CONCLUSION

The ALEXYS Clinical Analyzer in combination with a commercially available kit provides a standardised method for fast & reliable analysis of serotonin in urine.

## References

1. Smythe GA, Serotonin. In: Gray CH, v. James HT (Editors), *Hormones in Blood 3.*, Academic Press, New York, 125-134 (1983)
2. Cooper JR, Bloom RH, *The biochemical basis of neuropharmacology*, 5th Edition New York, Oxford University Press, (1986)
3. Stoward PJ, The possible role of 5-Hydroxytryptamine in Duchenne's muscular dystrophy. In: de Clerck F, Vanhoutte PM (Editors), *5-Hydroxytryptamine in peripheral reactions*, Raven Press, New York (1982)
4. Keller R, Greiling H, Gressner AM, *Lehrbuch der Klinischen Chemie und Pathobiochemie*, Schattauer Verlag, Stuttgart (1987)
5. Stahl SM, *Am. J. Psychiatry*, **26**, (1983) 140
6. Grahame-Smith DG, The carcinoid syndrome. In: Truelove SC, Lee E (Editors), *Topics in gastroenterology*, Blackwells, London (1972)
7. Kuhn DM, Lovenberg W, *Methods in biogenic amine research*, Elsevier Science B.V., (1983) Chapter 23
8. Manz B, *J. Clin. Biochem.*, **23**, (1985) 657
9. Mailman R, Kilts CD, *Clin. Chem.*, **31/11**, (1985) 1849-1854
10. Picard M, Olichon D, Gombert J, *J. Chromatogr.*, **341**, (1985) 445-451
11. Jouve J, Martineau J, Mariotte N, Barthelemy C, Muh JP, Lelord G, *J. Chromatogr.*, **378**, (1986) 437-443
12. Tagari PC, Boullin DJ, Davies CL, *Clin. Chem.* **30/1**, (1984) 131-135
13. R.T. Peaston and C. Weinkove, *Ann. Clin. Biochem.*, **41**, (2004) 17-38.
14. Recipe, Instruction manual for the serotonin in urine kit, version 3.4 (2005)

## PART NUMBERS AND CONFIGURATIONS

<b>180.0039C</b>	ALEXYS Clinical Analyzer
<b>110.4105</b>	VT03 3mm GC, salt bridge
<b>RE.7000</b>	ClinRep® complete kit , Serotonin in urine (for 100
<b>RE.7030</b>	ClinRep® Analytical column
<b>RE.8022</b>	ClinChek® urine controls, level I, II
<b>RE.7000</b>	ClinRep® complete kit , Serotonin in urine (for 100

Antec Leyden is an **ISO 9001:2008** certified company. For research purposes only. Specifications mentioned in this application note are subject to change without further notice. The information shown in this communication is solely to demonstrate the applicability of the ALEXYS for the analysis of metanephrines. The actual performance may be affected by factors beyond Antec Leyden's control.