

Amino acids using precolumn derivatization with DABS-Cl

A small amount of amino acids of protein and peptide can be quantitatively analyzed by vapor-phase hydrolyzing with hydrochloric acid, derivatizing with DABS-Cl (Dimethylamino azobenzene-sulfonyl chloride) and separated on reversed phase HPLC. DAB-Label is a kit to facilitate the troublesome work for derivatization. The obtained DAB-Labeled amino acids provide good stability and can be stored for quite a while.

Fig.1 shows the outline of the amino acid analysis by using DAB-Label, Fig. 2 shows the reaction formula of DABS-Cl and amino acid and Fig.3 shows the chromatogram of the derivatized representative amino acid samples.

Keywords: 1.Amino acid, 2.STD mixture, 3.ODS, 4.Vis, 5.DABS-Cl

Conditions:

| | |
|-------------------|--|
| Column: | Dabsylpak |
| Precolumn: | Dabsylpak-P |
| Eluent A: | 50mM NaH ₂ PO ₄ (pH6.6)/DMF(96/4) |
| Eluent B: | CH ₃ CN |
| Time(min) | 0.0 9.0 14.0 16.0 18.0 22.0 23.0 24.0 25.0 |
| A(%) | 85 70 59 49 46 36 10 10 85 |
| B(%) | 15 30 41 51 54 64 90 90 15 |
| 1Cycle | 40.0min |
| Wave length: | 465nm |
| Flow rate | 1.0ml/min |
| Column Temp: | 40 degree celsius |
| Sample: | Amino Acid Std. Soln. Type H (Wako Pure Chem) 50pmol/ul Histamine 500pmol/ul |
| Injection volume: | 20ul |




Fig.1 Outline of amino acid analysis by using DAB-Label




Fig. 2 Reaction formula of DABS-Cl and amino acids




Fig.3 Chromatogram of the derivatized representative amino acid samples