

## Analysis of Additives in Soft Drink Using Extreme Liquid Chromatography ( $\lambda$ -LC<sup>®</sup>)

### Introduction

Common additives in soft drink include caffeine, aspartame, benzoic acid, and acesulfame-K. Quality control demands to determine these additives. HPLC has been used for this analysis. An increase in the number of samples to be measured requires high-throughput analysis. In order to meet this requirement, we have applied extreme high pressure liquid chromatography system

( $\lambda$ -LC) equipped with a 2  $\mu$ m packing column to the analysis of these additives using high-pressure gradient method.

### Experimental

The system used for the measurement consists of two 3185PU pumps, 3080DG degasser, 3080MX mixing unit, 3067CO column oven 3070UV UV/Vis detector, 3059AS autosampler and a chromatography data system. Soft drink is degassed and filtered with 0.45  $\mu$ m membrane filter. A portion of 1  $\mu$ L is injected.

### Results and Discussion

Figure 1 shows the chromatogram of standard mixture of additives. Each component is sufficiently separated. The analysis time is within 2 minutes.  $\lambda$ -LC enables to shorten the analysis time significantly, leading to high-throughput analysis and to reduce consumption of solvent.




Figure 1 Chromatogram of the additives in soft drink

Components: 1=caffeine, 2=aspartame, 3=benzoic acid, 4=acesulfame K Injection volume: 1  $\mu$ L (filtrate) Column: X-Presspak C18S (2,1 mm I.D. x 50 mm, 2  $\mu$ m particles in diameter) column temperature: 40<sup>o</sup> C Mobile phase: A=50mM NaH<sub>2</sub>PO<sub>4</sub>(pH3.6), B=CH<sub>3</sub>OH Gradient condition: 0 min=A/B(90/10) > 1.5 min=A/B(50/50) > 2.0 min=A/B(50/50) > 2.05 min=A/B(90/10) Flow rate: 0.6 mL/min. Detection wavelength: 214 nm Pressure: 32 MPa-41MPa