

Analysis of alkylphenol in environmental water-1

Alkyl phenols are considered to be extrinsic endocrine disrupting chemicals that can disrupt the body's endocrine system ("environmental hormones"). In this report, tiny amounts of several alkyl phenols in environmental water were concentrated online using column switching and analyzed by reversed phase chromatography. Fig. 1. Shows the chromatograms of river water and river water spiked with 1 ppb of a standard sample.

Keywords: 1. alkyl phenol, 2. STD, river water, 3. ODS, 4. FL, 5. endocrine disrupting chemicals

STD

- | | |
|--------------------------|-----------------------|
| 1. 4,4-Dihydroxybiphenyl | 5. 4-Propylphenol |
| 2. 4-Ethylphenol | 6. 3-tert-Butylphenol |
| 3. Bisphenol A | 7. 4-tert-Butylphenol |
| 4. 4-Butoxyphenol | 8. 2-tert-Butylphenol |

Conditions:

Pretreatment Column: CrestPak C18T-5P
(4.6mm I.D. x 50mmL)

Pretreatment Eluent: H₂O

Switching time: 6 to 38min

Separation Column: CrestPak C18T-5
(4.6mm I.D. x 250mmL)

Eluent: A=CH₃OH / H₂O (60/40)

B=CH₃OH

Time(min)	A(%)	B(%)
0.0	100	0
35.0	100	0
35.1	0	100
50.0	0	100
50.1	100	0

1 cycle 70min

Flow rate: 1.0mL/min

Column temperature: 40 degree celsius

Sample: river water + STD(1ppb)

river water

Injection volume: 5mL

FL time program:	Time (min)	Ex(nm)	Em(nm)	Gain
	0.0	276	362	x1000
	13.0	280	310	x1000
	16.0	282	306	x1000
	18.2	290	328	x1000
	22.0	285	320	x1000
	27.0	298	302	x1000




Fig. 1. Chromatograms of river water