

Extraction of Acidic, Neutral and Basic Drugs from Urine using ISOLUTE® HCX Columns

Introduction

The following method addresses the extraction of all three classes of drugs of abuse from urine using non-polar and cation exchange retention mechanisms on the ISOLUTE HCX mixed mode column. This protocol was optimized to eliminate the use of dichloromethane in the elution of the acidic / neutral drug fraction.

Extraction Procedure

ISOLUTE® SPE Column: ISOLUTE® HCX 130 mg/10 mL, Part number 902-0013-H

Matrix: The matrix is aqueous, with relatively high ionic strength. There are many possible interference compounds, so this method uses a mixed retention mechanism, to allow extensive clean-up steps.

Pre-treatment: To urine (5 mL) in a clean dry test tube, add 0.1M phosphate buffer, pH 6.0 (2 mL). Vortex the sample, and check that the pH is 5.0–7.0. Adjust if necessary.

Condition: Condition the column with methanol (2 mL).

Equilibration: Rinse the column with deionized water (2 mL) followed by 0.1M phosphate buffer, pH 6.0 (3 mL).

Sample loading: Apply the sample to the column at a flow rate of 1–2 mL/min

Interference elution: Elute interferences with 0.1M phosphate buffer, pH 6.0 (1 mL) and dry the column under full vacuum for 5 mins. Add 1.0 M acetic acid (1 mL) and dry the column for a further 5 mins under full vacuum. Add hexane (1 mL).

Analyte elution: Acidic/neutral drug fraction. Elute the neutral and acidic drugs with ethyl acetate:hexane (25:75, v/v, 3 mL). Dry the column for 2 mins under full vacuum.

Interference elution. Elute interferences with methanol (6 mL) and dry the column under full vacuum for 2 mins. Basic drug elution. Elute the basic drugs with ethyl acetate: ammonium hydroxide (98:2, v/v, 3 mL).

Reagents:

- 1.0 M potassium hydroxide. Weigh potassium hydroxide (28 g) into a 500 mL volumetric flask. Add deionized water (400 mL), dissolve the potassium hydroxide and make up to the mark with deionized water.
- 0.1M phosphate buffer, pH 6.0. Weigh potassium hydrogen orthophosphate (13.61 g) into a one litre volumetric flask. Add deionized water (900 mL) and dissolve the potassium dihydrogen orthophosphate. Adjust the pH to 6.0 (+/- 0.1) with 1.0M potassium hydroxide, and make up to the mark with deionized water.
- 1.0M acetic acid. Add glacial acetic acid (5.75 mL) to a 100 mL volumetric flask. Add deionized water (70 mL), mix thoroughly, and make up to the mark with deionized water.
- 25:75 (v/v) ethyl acetate:hexane. Add ethyl acetate (25 mL) and hexane (75 mL) to a reagent bottle and mix thoroughly.
- 98:2 (v/v) ethyl acetate:ammonium hydroxide. Add concentrated ammonium hydroxide (2 mL) to a 100 mL volumetric flask. Add ethyl acetate (80 mL). Mix thoroughly, and make up to the mark with ethyl acetate.

Additional Information:

1. For reproducible and high recoveries, it is important to control the flow to 1–2 mL/min
2. The drying steps employed during the interference elution stage are important in order to maximize recoveries.
3. The ammonium hydroxide/ethyl acetate solution should be made up daily to obtain the required recoveries.
4. Previous # IST2006

NORTH AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
press (4) at the auto attendant
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
1-pointsupport@biotage.com

EUROPE

Main Office: +46 18 56 5900
Fax: +46 18 59 1922
Order Tel: +46 18 56 57 10
Order Fax: +46 18 56 57 05
order@eu.biotage.com

Japan

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com

Distributors

Please visit our Web site at
www.biotage.com
for contact details.