




# PROTOCOL MANAGEMENT SOFTWARE

MANAGE ALL INSTRUMENTS AND DISSOLUTION PROTOCOLS IN ONE DESKTOP VIEW

The screenshot displays the Protocol Manager software interface. On the left, there is a list of protocols under the heading "All Protocols", including "Acetohexamide Tablets", "Aspirin", "Auto-Start", "Autofilter Every Other", "Autofilter Every Sample", "Autofilter Every Third", "AutoMag Height Change", "Autowash Demo", "Betamethasone", "Caffeine", "Carbamazepine Tablets", "Demo - ADD and Autostart", "Demo - Speed Change", "Demo - Temp Change", "Dimenhydramine", "Dyphylline", "Dyphylline Tablets", "Enalapril Maleate", "Ethacrynic Acid Tablets", "FAST", "Ferrous Fumarate Tablets", "Flow Rate Change", and "Hydrochlorothiazide". A search bar and buttons for "Edit", "Delete", and "Print" are visible. The central area shows a search bar and a grid of instrument cards for "QC1", "QC2", and "R&D1". The "R&D1" section highlights the "V8 ENGINEERING" instrument. The right-hand panel provides detailed information for the selected "V8 ENGINEERING" instrument, including its serial number (0308-EB24), IP address (192.168.1.92), and a list of recent events and errors.

## NETWORKED PRODUCTIVITY

- Centralized management of users, groups, protocols, instruments, and security settings
- Automatically detect all instrument components and performance characteristics.
- Quickly drag-and-drop users, protocols, and methods to specific instruments.
- Convenient wizards help build complicated protocols; any conflicts are rapidly resolved through illustrated, context-sensitive help.
- Scalable from single local workstation to enterprise-wide across multiple lab sites



## ORDERING INFO

Protocol Manager software installation and operation qualification (IQ/OQ) services are also available. For further information visit [www.teledynehanson.com](http://www.teledynehanson.com) and contact your Teledyne Hanson representative directly or email us at [hansonsales@teledyne.com](mailto:hansonsales@teledyne.com).



9810 Variel Avenue  
Chatsworth, CA 91311 USA

+1 818.882.7266  
[www.teledynehanson.com](http://www.teledynehanson.com)

© Teledyne Hanson, a division of Teledyne Instruments, Inc. | 99-310-065 Rev. 07 Dec 2021