MSV-5800

UV/Vis/NIR Microscopic Spectrophotometer



Specifications

[Hardware]

Optical system: Single monochromator

Czerny-Turner mount Double beam type

Light source: Deuterium lamp, Halogen lamp

Option for light source: 150 W Xenon lamp

Wavelength range: 200 to 1600 nm

Wavelength accuracy: $\pm 0.3 \text{ nm} (656.1 \text{ nm})$

±1.5 nm (1312.2 nm)

Spectral bandwidth: 0.1, 0.2, 0.5, 1, 2, 5, 10 nm (UV/Vis region)

L2, L5, L10 nm (low stray light mode, UV/Vis region)

0.2, 0.4, 1, 2, 4, 10, 20 nm (NIR region)

L4, L10, L20 nm (low stray light mode, NIR region)

Scanning mode: Continuous scanning or step scanning

Detector: Photomultiplier tube, InGaAs photodiode

Mevasurement mode: Transmittance/reflectance

Sample observation: High-resolution built-in CMOS camera (3 million pixels), optical zoom, ATOS

(Aperture Through Optical System), LED illumination

Options for observation: Binocular, polarization observation unit, objective lens (10x, 20x)

Objective Mirrors: Cassegrain objective mirror (10x, 16x, 32x)

Automatic 4-position motor-driven objective revolver switching

Select one from 3 types of cassegrain mirror*1

Condensing Mirrors: Cassegrain condensing mirror (10x, 16x, 32x)

Manual replacement

With automatic correction of condensing mirror position

Select one from 3 types of cassegrain mirror*1

Aperture: Φ10, 20, 30, 50, 100, 200 μm

10 x 31.5, 10 x 50, 10 x 100, 31.5 x 10, 50 x 10, 100 x 10 μm

(when 16x Cassegrain objective mirror is used)

Sample stage: Manual stage

Moving distance: X-axis: 75 mm / Y-axis: 50 mm / Z-axis: 20 mm^{*2}

Options for sample stage: Automatic XYZ Stage, joystick

Moving distance: X-axis: 72 mm / Y-axis: 52 mm / Z-axis: 25 mm²

1 μm interval each for X, Y, Z-axis

Polarizer: Glan-Taylor prism

Automatic insertion/extraction to the light path, and automatic angular

setting

Analyzer (option): Glan-Taylor prism

Automatic insertion/extraction to the light path, and automatic angular

setting

Control panel: Objective mirror (lens) switching, measurement mode switching

(transmittance/reflectance), aperture switching, measurement start/stop, autofocus, auto-correction of transmittance focus, optical zoom, brightness control of observation light, sample compartment illumination (ON/OFF),

ATOS Illumination (ON/OFF)

Dimensions and weight: $740 \text{ (W)} \times 745 \text{ (D)} \times 630 \text{ (H)} \text{ mm, approx. } 111 \text{ kg}$

Power requirements: AC100 to 240 V, 50/60Hz, 340 VA

[Data processing]

Software: JASCO Spectra Manager Ver. 2

Operating system: Windows 10 Pro (64-bit)

Standard program

Measurement programs: [Microscopic Spectra Measurement] program, [Microscopic Fixed

Wavelength Measurement] program, [Microscopic Time Course Measurement] program (single-point measurement, image acquision)

Micro imaging analysis: Peak height, peak height ratio, peak area, peak area ratio, peak shift, full

width at half maximum

Spectra analysis: Film thickness analysis, color analysis, arithmetic, X and Y unit conversion,

derivatives, peak find, peak height/area, peak height/area ratio, full width

at half maximum, Kramers-Kronig transformation, FFT filter

Other programs: [Validation] program (jigs are required), [JASCO Canvas] program (print

layout tool), [Administrative Tools] program

Automatic XYZ stage operation:

(when mounting automatic XYZ stage)

Stage movement, lattice measurement, line measurement, multiplepoint measurement, fixed wavelength mapping measurement, autofocus,

multiple image acquision

Image display (mapping data): Color 3-D graph, bird's-eye view graph, contour map, color-coded map,

3-D spectrum graph, 2-D cross-sectional graph, RGB view, overlaying

view

Optional program: [Multi Layer Analysis] program, [Color Diagnosis Analysis] program

*1. MSV-5500/5700/5800 includes cassegrain objective mirror and cassegrain condensing mirror, whose magnifications are same.

*2. The moving distance of the Z-axis is limited by the magnification of the objective/condensing mirror or lens, or sample thickness.

The contents of this brochure are for reference and illustrative purposes only. Information, descriptions, and specifications in this publication are subject to change without notice. JASCO assumes no responsibility and will not be liable for any errors or omissions contained herein or for incidental, consequential damages or losses in connection with the furnishing, performance or use of this material.



